

# CK-7 Polyclonal Antibody

Catalog Number:D-AB-10469L

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


## Description

<b>Reactivity</b>	Human
<b>Immunogen</b>	Recombinant Human CK-7 protein expressed by E.coli
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification</b>	Antigen Affinity Purification
<b>Conjugation</b>	Unconjugated
<b>Formulation</b>	PBS with 0.02% sodium azide,1% protective protein and 50% glycerol,pH7.4


## Applications Recommended Dilution

<b>IHC</b>	1:500-1:1000
<b>IF</b>	1:50-1:200


## Data




Immunohistochemistry of paraffinembedded Human lung adenocarcinoma using CK-7 Polyclonal antibody at dilution of 1:1000.



Immunohistochemistry of paraffinembedded Human colon using CK-7 Polyclonal antibody at dilution of 1:500.



Immunohistochemistry of paraffinembedded Human kidney using CK-7 Polyclonal antibody at dilution of 1:500.



Immunofluorescence analysis of HeLa cells using CK-7 Polyclonal Antibody at dilution of 1:200

## Preparation & Storage

<b>Storage</b>	Store at -20°C Valid for 12 months. Avoid freeze / thaw cycles.
<b>Shipping</b>	The product is shipped with ice pack, upon receipt, store it immediately at the temperature recommended.

## Background

Keratins are a large family of proteins that form the intermediate filament cytoskeleton of epithelial cells, which are

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017

# CK-7 Polyclonal Antibody

Catalog Number:D-AB-10469L



classified into two major sequence types. Type I keratins are a group of acidic intermediate filament proteins,including K9–K23,and the hair keratins Ha1–Ha8. Type II keratins are the basic or neutral counterparts to the acidic type I keratins,including K1–K8,and the hair keratins,Hb1–Hb6. KRT7,also named as cytokeratin 7,is one member of type II basic cytokeratin. It is specifically expressed in the simple epithelia lining the cavities of the internal organs and in the gland ducts and blood vessels,and their neoplasms. KRT7 is marker of epithelial tissues,but not present in carcinomas of stratified squamous cell origin.This antibody is specifically against KRT7.

---

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Web: [www.elabscience.com](http://www.elabscience.com)

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)